Experiment 1: Determination of free CO₂ using N/50 NaCO₃

- 1) What is the Indicator used in this reaction? Ans: Phenolphthalein.
- 2) What is the End-point of this reaction? Ans: Colourless to **Pink**.
- 3) What is Normality. What is its S.I unit? Ans: The number of gram or mole equivalents of solute present in one litre of a solution. The SI unit for normality is moles per liter (mol/L). Its denoted by "N".
- 4) What is strength of a solution? Ans: The strength of the solution is defined as the amount of solute in gram dissolved per liter of the solution.

S.I unit: mg/L or ppm.

5) What type of titration is performed in this experiment?

Ans: Acid-Base titration.

6) Draw the structure of Phenolphthalein.

7) Is CO₂ Acidic or Basic? Ans: Acidic.

8) What do you mean by free CO₂ in water? Ans: "Free CO₂ in water" refers to the dissolved carbon dioxide gas present in a water sample that is not chemically bound or in equilibrium with other dissolved species. When carbon dioxide (CO₂) is exposed to water, it can dissolve and form carbonic acid (H2CO₃), which can further dissociate into bicarbonate ions (HCO₃-) and hydrogen ions (H+). The term "free CO₂" specifically refers to the dissolved CO₂ that remains in its gaseous form in the water.

Experiment 2: Determination of Moisture in Soil.

- 1) What is the use of Desiccator? Ans: A desiccator is a specialized piece of laboratory equipment used for the drying or preservation of moisture-sensitive substances. It consists of a sealed container made of glass or plastic with a removable lid and an airtight seal. The interior of the desiccator contains a drying agent such as silica gel or anhydrous CaCl₂.
- 2) What type of method is used for measuring the content of water in soil? Ans: Gravimetric.
- 3) What is the temperature to be set in the Oven during experiment?

Ans: 105° C.

Experiment 3: Determination of Ammonia NH_3 in water using N/50 H_2SO_4

1) What form ammonia is present in water? Ans: Ammonium Hydroxide (NH₄OH).

2) What type of titration is performed during the experiment?

Ans: Acid-Base.

3) What is the indicator used in this reaction? Ans: Methyl-Red.

4) What is the End-point of this reaction? Ans: Yellow to **Red**.

5) Draw the structure of methyl red. Ans:

Experiment 4: Determination TDS in water.

- 1) What is TDS (Total Dissolved Solids)? Ans: Dissolved solids refers to any minerals, salts, cation or anions dissolved in water, which varies between 20-1000 ppm.
- 2) What is the temperature to be set in oven for this experiment?

Ans: 180°-190° C.

Experiment 5: TRC in Water using N/50 Na₂S₂O₃ and Conc. HCL.

1) What kind of titration is performed in this experiment? Ans Idometric as I₂ is liberated.

2) How is residual chlorine estimated? Ans: Residual Chlorine is estimated by oxidation of KI with RC with liberation of I₂.

3) What is indicator used? Ans: Starch. Blue-Black.

- 4) Colours shown: Brown to straw yellow to blue-black to colourless.
- 5) Why we don't use pipette for this experiment. Ans: Because chlorine vapours are harmful so the solution is not sucked into pipette. Instead Measuring cylinder is used.
- 6) What is residual Chlorine? Ans: Residual chlorine refers to the amount of chlorine that remains in water or a solution after disinfection or chlorination processes have

occurred. Chlorine is widely used as a disinfectant in water treatment to kill or inactivate harmful microorganisms, including bacteria, viruses, and parasites.

Experiment 6: <u>Dissolved Oxygen in water using</u> <u>N/40 Na₂S₂O₃ and Conc. H₂SO₄</u>

1) What are the oxygen carrier in this experiment?

Ans: MnSO₂.

2) What kind of titration is performed in this experiment? Ans Idometric as I₂ is liberated.

- 3) Why we check oxygen levels for pollution? Ans: Oxygen itself is not a pollution, but its absence is an indicator of several types of pollution.
- 4) What is azide in Alkaline Iodide Azide? Ans: N_3^- .
- 5) What is used remove nitrites (NO₂⁻)? Ans: Sodium Azide.

6) Draw structure of $C_6H_{10}O_5$ (Starch)? What is its colour?

Ans: Blue-Black colour.

7) Why we didn't add starch indicator at the beginning of the reaction?

Ans: Because it will reach will other reagents.

- 8) In what state oxygen is present in water? Ans: Oxygen is present in its molecular form. So iodine is not capable of reacting.
- 9) What is the difference between end point and equivalence point?
- 10) What is Hypo solution?

Ans: Na₂S₂O₃.

Experiment 7 and 8: Determination of PH and Conductivity of water.

1) What is the S.I unit of conductivity? Ans: Siemens per meter.